




WES and WGS Planning and Best Practices



General considerations to increase the likelihood of high quality data



Technology comparison

	Whole-exome-seq (WES)	Whole-genome-seq (WGS)
Intended capture	Exons (~2% of genome) Requires an exon capture kit	Whole-genome
Types of calls	Germline and somatic (~5% detection limit of somatic)	Germline and somatic (~10-15% detection limit of somatic)
Typical coverage	100-300X, depth is highly variable	30X, typically fairly uniform
Intended variants	SNPs and indels within exons	CNA, gene fusions, SNPs and indels in high frequency
Cost	Lower	More costly than WES because of sequencing needed

Alternative technologies

Duplex-seq

Tag-based error correction method to improve sequencing accuracy of select target regions
Capable of detecting mutations <0.1%

Genotyping arrays or low-pass WGS


Lower-cost alternative if interested in only germline alterations
Does not cover every base

Custom gene panels

Useful if only specific regions are of interest

Long-read sequencing

Provides better resolution of more challenging genomic regions and structural variants



Quality control

Material evaluation

Assess concentration, fragment size, DIN

Library evaluation

Check size and quality (e.g. bioanalyzer, tapestation)
Pilot studies and shallow sequencing can identify issues

Data evaluation

Check QC metrics (e.g. depth, contamination, coverage regions)
Validate important variants

Decide on sequencing

Evaluate material and library QC and seek expert consultation
Go/no-go decision



Contact the Biostatistics & Bioinformatics Shared Resource (BBSR) for cancer-related analysis

James.Costello@CUAnschutz.edu, Bioinformatics@CUAnschutz.edu



Contact CUBE for all non-cancer-related projects

