




Spatial Transcriptomics Planning and Best Practices



General considerations to increase the likelihood of high quality data



Spatial transcriptomics: RNA mapped to tissue coordinates

Spatial Considerations

What are your **spatial** questions?

These platforms are most applicable for the study of gene expression and how it relates to tissue architecture, cell niches, cell-cell interactions, and microenvironments.

Gene count depth and breadth is generally lower than scRNAseq or similar methods. This may affect the ability to reliably detect genes with low expression.

Sample Handling

Most platforms can handle FFPE and fresh frozen samples.

For best results use tissue stored for < 1-2 years.

Confirm tissue RNA quality using RNA integrity (RIN) or RNAscope >80% positivity with positive control probes.



Design Parameters

Always mix control and treatment samples on slides and single runs to allow for instrument level normalization of signal detection. Biological replication is expensive but always encouraged. Direct comparisons between separate runs and longitudinal studies are not currently recommended due to long instrument processing times and reagent lot differences.

Analysis / QC

Review staining, image quality and segmentation outputs to confirm accurate cell boundaries. Some exploration is possible immediately after run (AtoMx, Loupe Browser, Xenium Explorer). Most questions will require experience with a coding environment (R/Python) and specialized packages (Giotto, SpatialExperiment, Scanpy/Squidpy). This form of analysis is lengthy and requires iterative collaboration between analyst and investigator.

Consideration of platform choice

Assay Type

FFPE and Fresh Frozen samples are accepted

Assay readout types include:

- imaging-based targeted *in situ* RNA with fluorescent probes
- oligo probe based sequencing
- 3' polyA capture sequencing

Pre-designed fluorescent probe panels range from ~50 to 5000 genes or whole transcriptome.

Detection

How transcripts are associated with biological features (i.e. cell segmentation) can differ based on spatial method:

- grid based coordinates
- a user defined region of interest
- membrane staining
- nuclear distance

The grid overlay strategy is 2µm, 8µm, and 16µm resolutions.

Transcript localization can be ~50 to 120nm precision.

Notes

Custom probe additions are typically recommended.

Some platforms offer additional post-hoc staining and protein detection as addons.

Sequencing based platforms can be paired with VDJ and long read sequencing.

Imaging platforms have long slide processing times (~2 to 6 days or more).

Downstream analysis outside of basic exploration will require significant computation, large data storage capabilities, and analyst support.

Platforms

Consider the information on this flyer and then meet with an analyst and the appropriate campus core to discuss the solution that best fits your needs.

CU Anschutz instruments:

Genomics Shared Resource:

- genomics.core@cuanschutz.edu
- 10X Genomics Visium HD
- Bruker/Nanostring GeoMx
- Bruker/Nanostring CosMx

Human Immune Monitoring Shared Resource (HIMSR):

- Angela.Minic@cuanschutz.edu
- 10X Genomics Xenium

Spatial Transcriptomics Office Hours; Hosted jointly by HIMSR and GSR cores: Every Friday from 10:00 am - 11:00 am MST. <https://ucdenver.zoom.us/j/97848764350>

Contact the Biostatistics & Bioinformatics Shared Resource (BBSR) for cancer-related analysis

James.Costello@CUAnschutz.edu, Bioinformatics@CUAnschutz.edu



Contact CUBE for all non-cancer-related projects

